



Recombinant HCV NS3 [His] (DAG4635)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	> 98% as determined by SDS-PAGE
Conjugate	His
Applications	ELISA
Format	Liquid
Size	1 mg
Buffer	25 mM Tris-HCl, 0.1% SDS, pH 8.3
Preservative	None
Storage	Store at -20°C

BACKGROUND

Introduction	The hepatitis C virus (HCV) is a small (55–65 nm in size), enveloped, positive-sense single-stranded RNA virus of the family Flaviviridae. The hepatitis C virus is the cause of hepatitis C and some cancers such as liver cancer (hepatocellular carcinoma, abbreviated HCC) and lymphomas in humans.
Keywords	Hepatitis C virus; HCV; NS3