



Native Aspergillus Giganteus Alpha-Sarcin (DAG4644)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Aspergillus giganteus Alpha-Sarcin
Species	Aspergillus
Conjugate	Unconjugated
Applications	immunogen
Format	freeze dried
Concentration	trace amount of the dimeric form.
Preservative	None
Storage	2-8°C short term, -20°C long term

BACKGROUND

Introduction	Aspergillus giganteus Alpha-Sarcin
Keywords	Alpha-sarcin; a-sarcin; Aspergillus giganteus Alpha-Sarcin; Aspergillus; Aspergillus giganteus