



Human C8 Protein (DAG4673)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Human C8 Protein
Species	Human
Purity	> 90% by SDS-PAGE
Conjugate	Unconjugated
Applications	immunogen
Format	Liquid
Concentration	1.0 mg/ml (see Certificate of Analysis for actual concentration)
Size	0.25 mg
Buffer	10 mM sodium phosphate, 145 mM NaCl, pH 7.3
Preservative	None
Storage	2-8°C short term, -20°C long term

BACKGROUND

Keywords	C8
----------	----