



# Tobramycin [HSA-Biotin] (DAG469S)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Antigen Description	Tobramycin, covalently linked to biotinylated human serum albumin
Species	N/A
Purity	chromatography, dialysis
Conjugate	HSA-Biotin
Applications	e.g. allergy test, immunoassays etc.
Molecular Weight	app. 65 kDa
Format	Liquid
Concentration	1 mg/mL, conjugation ratio is app. 1:10 - 1:40.
Size	1 mg
Buffer	PBS
Preservative	0.05% sodium azide
Storage	Store at 4°C.

## BACKGROUND

Introduction	Tobramycin is an aminoglycoside antibiotic used to treat infections caused by many different bacteria. Tobramycin is usually administered by intravenous infusion, intramuscular injection, or inhalation.
Keywords	drug-HSA-Biotin conjugate, drug-conjugate, tobramycin HSA-Biotin conjugate