



Recombinant CMV (gB) [GST] (DAG4719)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	CMV-gB (residues 11-67), GST-tagged
Specificity	Immunoreactive with sera of CMV-infected individuals.
Species	CMV
Purity	> 95% by SDS-PAGE
Conjugate	GST
Applications	Antigen in ELISA and Western blots, excellent antigen for detection of CMV with minimal specificity problems.
Format	Liquid. Supplied in 25mM Tris-HCl pH 8.0, 1mM EDTA and 50% glycerol.
Preservative	None
Storage	2-8°C short term, -20°C long term

BACKGROUND

Introduction	Cytomegalovirus (from the Greek cyto-, "cell", and -megalo-, "large") is a viral genus of the viral family known as Herpesviridae or herpesviruses. It is typically abbreviated as CMV.
Keywords	Cytomegalovirus gB; CMV gB; Cytomegalovirus; CMV