



# Recombinant CMV Pp38 [GST] (DAG4731)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Product Overview</b>	Recombinant Cytomegalovirus pp38 (UL80a), GST-tagged
<b>Specificity</b>	Immunoreactive with sera of CMV-infected individuals.
<b>Species</b>	CMV
<b>Purity</b>	> 95%, as determined by SDS-PAGE.
<b>Conjugate</b>	GST
<b>Applications</b>	CMV pp38 antigen is suitable for ELISA and Western blots, excellent antigen for detection of CMV with minimal specificity problems.
<b>Format</b>	Lyophilized from 25mM Tris-HCl pH 7.2, 1mM EDTA and 50% glycerol.
<b>Preservative</b>	None
<b>Storage</b>	2-8°C short term, -20°C long term

## BACKGROUND

<b>Introduction</b>	Cytomegalovirus (from the Greek cyto-, "cell", and -megalo-, "large") is a viral genus of the viral family known as Herpesviridae or herpesviruses. It is typically abbreviated as CMV.
<b>Keywords</b>	Cytomegalovirus Pp38; CMV Pp38; Cytomegalovirus; CMV