



Simvastatin [HSA-Biotin] (DAG480S)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Simvastatin, covalently linked to biotinylated human serum albumin
Species	N/A
Purity	chromatography, dialysis
Conjugate	HSA-Biotin
Applications	e.g. allergy test, immunoassays etc.
Molecular Weight	app. 65 kDa
Format	Liquid
Concentration	1 mg/mL, conjugation ratio is app. 1:10 - 1:40.
Size	1 mg
Buffer	PBS
Preservative	0.05% sodium azide
Storage	Store at 4°C.

BACKGROUND

Introduction	Simvastatin is a lipid-lowering agent derived synthetically from a fermentation product of the fungus <i>Aspergillus terreus</i> . Hydrolyzed in vivo to an active metabolite, simvastatin competitively inhibits hepatic hydroxymethyl-glutaryl coenzyme A (HMG-CoA) reductase, the enzyme which catalyzes the conversion of HMG-CoA to mevalonate, a key step in cholesterol synthesis. This agent lowers plasma cholesterol and lipoprotein levels, and modulates immune responses by suppressing MHC II (major histocompatibility complex II) on interferon gamma-stimulated,
--------------	---

antigen-presenting cells such as human vascular endothelial cells.

Keywords

drug-HSA-Biotin conjugate, drug-conjugate, simvastatin HSA-Biotin conjugate
