



Human Ribosomal Phosphoprotein P2 [His] (DAG4845)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Human
Purity	>80% as determined by SDS-PAGE
Conjugate	His
Applications	WB
Molecular Weight	13 kDa
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	50 µg
Preservative	None
Storage	2-8°C short term, -20°C long term

BACKGROUND

Introduction	The ribosomal phosphoproteins, also called P protein antigens, are associated with the large ribosomal subunit and therefore are antigenic targets with a cytoplasmic localization. Three P proteins have been described: P0 with a molecular weight of 35 kDa, P1 (19 kDa) and P2 (17 kDa). Autoantibodies against ribosomal P proteins are present in 10% of SLE patients. If anti-ribosomal P antibodies were to occur in the absence of other typical SLE associated autoantibodies, they may account for some patients with so-called ANA-negative lupus. It has been reported that lupus patients positive for anti-ribosomal P autoantibodies have a high
---------------------	---

frequency of CNS involvement, suggesting a marker use for these antibodies.

Keywords

Ribosomal phosphoprotein; P0; P1; P2
