



Human PM/ScI100 [His] (DAG4858)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Human
Purity	>70% as determined by SDS-PAGE
Conjugate	His
Applications	WB
Molecular Weight	102 kDa
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	50 µg
Preservative	None
Storage	2-8°C short term, -20°C long term

BACKGROUND

Introduction	PM/ScI 100 is a nucleolar protein and represents one crucial component of the so-called "PM/ScI complex" or "nuclear exosome", which is an RNA-processing complex consisting of a number of exoribonucleases (not to be confused with certain multivesicular bodies which are shed by cultured cells into the culture fluid and are also called exosomes)
Keywords	Polymyositis/scleroderma; PM/ScI