



Human M2 [His] (DAG4875)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|---------------|--|
| Species | Human |
| Purity | >75% as determined by SDS-PAGE |
| Conjugate | His |
| Applications | WB |
| Format | Liquid |
| Concentration | Batch dependent - please inquire should you have specific requirements |
| Size | 50 µg |
| Preservative | None |
| Storage | 2-8°C short term, -20°C long term |

BACKGROUND

| | |
|--------------|--|
| Introduction | Mitochondrial 2-oxo acid dehydrogenase multienzyme complexes play a pivotal role in the regulation of the cellular redox state |
| Keywords | 2 oxo acid dehydrogenase complex; M2 |