



Piroxycam [HSA-Biotin] (DAG488S)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|----------------------------|--|
| Antigen Description | Piroxycam, covalently linked to biotinylated serum albumin |
| Species | N/A |
| Purity | chromatography, dialysis |
| Conjugate | HSA-Biotin |
| Applications | e.g. allergy test, immunoassays etc. |
| Molecular Weight | app. 65 kDa |
| Format | Liquid |
| Concentration | 1 mg/mL, conjugation ratio is app. 1:10 - 1:40. |
| Size | 1 mg |
| Buffer | PBS |
| Preservative | 0.05% sodium azide |
| Storage | Store at 4°C. |

BACKGROUND

| | |
|---------------------|--|
| Introduction | Piroxicam is a nonsteroidal anti-inflammatory drug of the oxicam class used to relieve the symptoms of painful inflammatory conditions like arthritis. |
| Keywords | drug-HSA-Biotin conjugate, drug-conjugate, piroxycam HSA-Biotin conjugate |