



# Nystatine [HSA-Biotin] (DAG499S)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Antigen Description	Nystatine, covalently linked to biotinylated human serum albumin
Species	N/A
Purity	chromatography, dialysis
Conjugate	HSA-Biotin
Applications	e.g. allergy test, immunoassays etc.
Molecular Weight	app. 65 kDa
Format	Liquid
Concentration	1 mg/mL, conjugation ratio is app. 1:10 - 1:40.
Size	1 mg
Buffer	PBS
Preservative	0.05% sodium azide
Storage	Store at 4°C.

## BACKGROUND

Introduction	Nystatin is a topical and oral antifungal agent with activity against many species of yeast and candida albicans, which is used largely to treat skin and oropharyngeal candidiasis. Nystatin is not absorbed orally and has not been linked to drug induced liver injury.
Keywords	drug-HSA-Biotin conjugate, drug-conjugate, Nystatine HSA-Biotin conjugate