



Neomycin [HSA-Biotin] (DAG500S)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Neomycin, covalently linked to biotinylated human serum albumin
Species	N/A
Purity	chromatography, dialysis
Conjugate	HSA-Biotin
Applications	e.g. allergy test, immunoassays etc.
Molecular Weight	app. 65 kDa
Format	Liquid
Concentration	1 mg/mL, conjugation ratio is app. 1:10 - 1:40.
Size	1 mg
Buffer	PBS
Preservative	0.05% sodium azide
Storage	Store at 4°C.

BACKGROUND

Introduction	Neomycin is an aminoglycoside broad spectrum antibiotic effective against both gram positive and gram negative bacteria. It interferes with protein synthesis in sensitive bacterial cells such as species of Proteus and Staphylococcus.
Keywords	drug-HSA-Biotin conjugate, drug-conjugate, Neomycin HSA-Biotin conjugate