



n-butyl-p-aminobenzoate [HSA-Biotin] (DAG503S)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	n-butyl-p-aminobenzoate, covalently linked to biotinylated human serum albumin
Species	N/A
Purity	chromatography, dialysis
Conjugate	HSA-Biotin
Applications	e.g. allergy test, immunoassays etc.
Molecular Weight	app. 65 kDa
Format	Liquid
Concentration	1 mg/mL, conjugation ratio is app. 1:10 - 1:40.
Size	1 mg
Buffer	PBS
Preservative	0.05% sodium azide
Storage	Store at 4°C.

BACKGROUND

Introduction	N-BUTYL-P-AMINO BENZOATE is an aminophenyl ester derivative. Amines are chemical bases. They neutralize acids to form salts plus water.
Keywords	drug-HSA-Biotin conjugate, drug-conjugate, n-butyl-p-aminobenzoate HSA-Biotin conjugate