



Metoprolol [HSA-Biotin] (DAG506S)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Metoprolol, covalently linked to biotinylated human serum albumin
Species	N/A
Purity	chromatography, dialysis
Conjugate	HSA-Biotin
Applications	e.g. allergy test, immunoassays etc.
Molecular Weight	app. 65 kDa
Format	Liquid
Concentration	1 mg/mL, conjugation ratio is app. 1:10 - 1:40.
Size	1 mg
Buffer	PBS
Preservative	0.05% sodium azide
Storage	Store at 4°C.

BACKGROUND

Introduction	Metoprolol is a cardioselective competitive beta-1 adrenergic receptor antagonist with antihypertensive properties and devoid of intrinsic sympathomimetic activity. Metoprolol antagonizes beta 1-adrenergic receptors in the myocardium, thereby reducing the rate and force of myocardial contraction leading to a reduction in cardiac output. This agent may also reduce the secretion of renin with subsequent reduction in levels of angiotensin II thereby preventing vasoconstriction and aldosterone secretion.
--------------	---

Keywords

drug-HSA-Biotin conjugate, drug-conjugate, Metoprolol HSA-Biotin conjugate
