



# Recombinant HIV type 1 p24 (Antigen) (DAG515)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Purity	> 90% as determined by SDS-PAGE
Conjugate	TBD
Applications	ELISA, LFIA, WB, CLIA
Format	Liquid
Size	1 mg
Buffer	20mM PBS, 7.4
Storage	Store at -20°C

## BACKGROUND

Keywords	HIV1; p24; Human immunodeficiency virus
----------	---