



Cotrimoxazol [HSA-Biotin] (DAG528S)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Cotrimoxazol, covalently linked to biotinylated human serum albumin
Species	N/A
Purity	chromatography, dialysis
Conjugate	HSA-Biotin
Applications	e.g. allergy test, immunoassays etc.
Molecular Weight	app. 65 kDa
Format	Liquid
Concentration	1 mg/mL, conjugation ratio is app. 1:10 - 1:40.
Size	1 mg
Buffer	PBS
Preservative	0.05% sodium azide
Storage	Store at 4°C.

BACKGROUND

Introduction	Co-trimoxazole is used to treat certain bacterial infections, such as pneumonia (a lung infection), bronchitis (infection of the tubes leading to the lungs), and infections of the urinary tract, ears, and intestines. It also is used to treat 'travelers' diarrhea. Co-trimoxazole is a combination of trimethoprim and sulfamethoxazole and is in a class of medications called sulfonamides. It works by stopping the growth of bacteria. Antibiotics will not kill viruses that can cause colds, flu, or other viral infections.
--------------	---

Keywords

drug-HSA-Biotin conjugate, drug-conjugate, Cotrimoxazol HSA-Biotin conjugate
