



Ampicillin [HSA-Biotin] (DAG551S)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Ampicillin, covalently linked to biotinylated human serum albumin
Species	N/A
Purity	chromatography, dialysis
Conjugate	HSA-Biotin
Applications	e.g. allergy test, immunoassays etc.
Molecular Weight	app. 65 kDa
Format	Liquid
Concentration	1 mg/mL, conjugation ratio is app. 1:10 - 1:40.
Size	1 mg
Buffer	PBS
Preservative	0.05% sodium azide
Storage	Store at 4°C.

BACKGROUND

Introduction	Ampicillin is a broad-spectrum, semi-synthetic, beta-lactam penicillin antibiotic with bactericidal activity. Ampicillin binds to and inactivates penicillin-binding proteins (PBP) located on the inner membrane of the bacterial cell wall. Inactivation of PBPs interferes with the cross-linkage of peptidoglycan chains necessary for bacterial cell wall strength and rigidity. This interrupts bacterial cell wall synthesis and results in the weakening of the bacterial cell wall and causes cell lysis. Ampicillin is stable against hydrolysis by a variety of beta-lactamases, therefore, can
--------------	--

be used in wide range of gram-positive and -negative infections.

Keywords

drug-HSA-Biotin conjugate, drug-conjugate, Ampicillin HSA-Biotin conjugate
