



Amoxicillin [HSA-Biotin] (DAG553S)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Amoxicillin, covalently linked to biotinylated human serum albumin
Species	N/A
Purity	chromatography, dialysis
Conjugate	HSA-Biotin
Applications	e.g. allergy test, immunoassays etc.
Molecular Weight	app. 65 kDa
Format	Liquid
Concentration	1 mg/mL, conjugation ratio is app. 1:10 - 1:40.
Size	1 mg
Buffer	PBS
Preservative	0.05% sodium azide
Storage	Store at 4°C.

BACKGROUND

Introduction	Amoxicillin Anhydrous is the anhydrous form of a broad-spectrum, semisynthetic aminopenicillin antibiotic with bactericidal activity. Amoxicillin binds to and inactivates penicillin-binding proteins (PBPs) located on the inner membrane of the bacterial cell wall. Inactivation of PBPs interferes with the cross-linkage of peptidoglycan chains necessary for bacterial cell wall strength and rigidity. This interrupts bacterial cell wall synthesis and results in the weakening of the bacterial cell wall and causes cell lysis.
--------------	--

Keywords

drug-HSA-Biotin conjugate, drug-conjugate, Amoxicillin HSA-Biotin conjugate
