



Allopurinol [HSA-Biotin] (DAG558S)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Allopurinol, covalently linked to biotinylated human serum albumin
Species	N/A
Purity	chromatography, dialysis
Conjugate	HSA-Biotin
Applications	e.g. allergy test, immunoassays etc.
Molecular Weight	app. 65 kDa
Format	Liquid
Concentration	1 mg/mL, conjugation ratio is app. 1:10 - 1:40.
Size	1 mg
Buffer	PBS
Preservative	0.05% sodium azide
Storage	Store at 4°C.

BACKGROUND

Introduction	Allopurinol is a structural isomer of hypoxanthine. Allopurinol inhibits xanthine oxidase, an enzyme that converts oxypurines to uric acid. By blocking the production of uric acid, this agent decreases serum and urine concentrations of uric acid, thereby providing protection against uric acid-mediated end organ damage in conditions associated with excessive production of uric acid, i.e. the massive cell lysis associated with the treatment of some malignancies.
--------------	--

Keywords

drug-HSA-Biotin conjugate, drug-conjugate, Allopurinol HSA-Biotin conjugate
