



## Alcuronium [HSA-Biotin] (DAG559S)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

Antigen Description	Alcuronium, covalently linked to biotinylated human serum albumin
Species	N/A
Purity	chromatography, dialysis
Conjugate	HSA-Biotin
Applications	e.g. allergy test, immunoassays etc.
Molecular Weight	app. 65 kDa
Format	Liquid
Concentration	1 mg/mL, conjugation ratio is app. 1:10 - 1:40.
Size	1 mg
Buffer	PBS
Preservative	0.05% sodium azide
Storage	Store at 4°C.

### BACKGROUND

Introduction	Alcuronium is an indole alkaloid of the curare family. A neuroblocker, it is often used in chloride form as an anesthesia adjuvant. It has a role as a muscle relaxant.
Keywords	drug-HSA-Biotin conjugate, drug-conjugate, Alcuronium HSA-Biotin conjugate