



## Carbadox [KLH] (DAG5613)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

<b>Antigen Description</b>	Carbadox is an antibacterial agent used in swine, which can promote growth of pigs. It also as a effective control of swine dysentery (vibronic dysentery, bloody scours, or hemorrhagic dysentery) as well as bacterial swine enteritis (salmonellosis or necrotic enteritis caused by Salmonella enterica). Besides, it aid in the prevention of migration and establishment of large roundworm (Ascaris suum) infections and establishment of nodular worm (Oesophagostomum) infectionsand. But it has been banned in Canada, Australia, and the European Union.
<b>Species</b>	N/A
<b>Conjugate</b>	KLH
<b>Applications</b>	The antigen can be used in ELISA assay.
<b>Format</b>	Liquid
<b>Concentration</b>	Batch dependent - please inquire should you have specific requirements.
<b>Size</b>	1 mg
<b>Buffer</b>	0.01M pH7.4 PBS
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.