



Praziquantel [BSA] (DAG5620)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description

Praziquantel is a pyrazinoisoquinoline derivative with anthelmintic property. It is an antihelmintic agent with activity against a broad spectrum of trematodes and cestodes that is used predominantly in the therapy of schistosomiasis, liver flukes, and cysticercosis. Praziquantel increases the permeability of the tegument of susceptible worms, resulting in an influx and increase in intra-tegumental calcium leading to rapid contractions and paralysis of the worm's musculature through a subsequent increase in levels of calcium in the sarcoplasmic reticulum. In addition, vacuolization of the tegumental syncytium and blebbing results in tegument disintegration, leads to antigen exposure and elicit host defense responses to the worm.

Species	N/A
Conjugate	BSA
Applications	The antigen can be used in ELISA assay.
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.