



Moroxydine hydrochloride [KLH] (DAG5628)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Moroxydine hydrochloride is an antiviral drug that was developed for the treatment of influenza virus and herpes virus infection. It has potential applications against a number of RNA and DNA viruses. Structurally moroxydine is a heterocyclic biguanidine.
Species	N/A
Conjugate	KLH
Applications	The antigen can be used in ELISA assay.
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.