



Imidocarb [KLH] (DAG5631)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Imidocarb is a urea derivative used in veterinary medicine as an antiprotozoal agent for the treatment of infection with Babesia (babesiosis) and other parasites.
Species	N/A
Conjugate	KLH
Applications	The antigen can be used in ELISA assay.
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.