



## Imidocarb [BSA] (DAG5632)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

|                            |                                                                                                                                                                    |
|----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Antigen Description</b> | Imidocarb is a urea derivative used in veterinary medicine as an antiprotozoal agent for the treatment of infection with Babesia (babesiosis) and other parasites. |
| <b>Species</b>             | N/A                                                                                                                                                                |
| <b>Conjugate</b>           | BSA                                                                                                                                                                |
| <b>Applications</b>        | The antigen can be used in ELISA assay.                                                                                                                            |
| <b>Format</b>              | Liquid                                                                                                                                                             |
| <b>Concentration</b>       | Batch dependent - please inquire should you have specific requirements.                                                                                            |
| <b>Size</b>                | 1 mg                                                                                                                                                               |
| <b>Buffer</b>              | 0.01M pH7.4 PBS                                                                                                                                                    |
| <b>Storage</b>             | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.                                                                              |