



## Triazophos [KLH] (DAG5634)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

<b>Antigen Description</b>	Triazophos is derives from 1-phenyl-1H-1,2,4-triazol-3-ol. It is an organic thiophosphate and an organothiophosphate insecticide, and plays a role as an EC 3.1.1.7 (acetylcholinesterase) inhibitor, an acaricide, an agrochemical and a nematicide.
<b>Species</b>	N/A
<b>Conjugate</b>	KLH
<b>Applications</b>	The antigen can be used in ELISA assay.
<b>Format</b>	Liquid
<b>Concentration</b>	Batch dependent - please inquire should you have specific requirements.
<b>Size</b>	1 mg
<b>Buffer</b>	0.01M pH7.4 PBS
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.