



Triazophos [BSA] (DAG5635)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Triazophos is derived from 1-phenyl-1H-1,2,4-triazol-3-ol. It is an organic thiophosphate and an organothiophosphate insecticide, and plays a role as an EC 3.1.1.7 (acetylcholinesterase) inhibitor, an acaricide, an agrochemical and a nematicide.
Species	N/A
Conjugate	BSA
Applications	The antigen can be used in ELISA assay.
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.