



Isoprocarb [OVA] (DAG5636)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Isoprocarb is a carbamate ester, derives from a methylcarbamic acid and a 2-isopropylphenol. It has a role as an EC 3.1.1.7 (acetylcholinesterase) inhibitor, a carbamate insecticide and an agrochemical used on rice, cocoa, sugar cane, vegetables and other crop.
----------------------------	---

Species	N/A
----------------	-----

Conjugate	OVA
------------------	-----

Applications	The antigen can be used in ELISA assay.
---------------------	---

Format	Liquid
---------------	--------

Concentration	Batch dependent - please inquire should you have specific requirements.
----------------------	---

Size	1 mg
-------------	------

Buffer	0.01M pH7.4 PBS
---------------	-----------------

Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
----------------	---