



Isoprocarb [KLH] (DAG5637)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Isoprocarb is a carbamate ester, derives from a methylcarbamic acid and a 2-isopropylphenol. It has a role as an EC 3.1.1.7 (acetylcholinesterase) inhibitor, a carbamate insecticide and an agrochemical used on rice, cocoa, sugar cane, vegetables and other crop.
Species	N/A
Conjugate	KLH
Applications	The antigen can be used in ELISA assay.
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.