



## Isoprocarb [BSA] (DAG5638)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

<b>Antigen Description</b>	Isoprocarb is a carbamate ester, derives from a methylcarbamic acid and a 2-isopropylphenol. It has a role as an EC 3.1.1.7 (acetylcholinesterase) inhibitor, a carbamate insecticide and an agrochemical used on rice, cocoa, sugar cane, vegetables and other crop.
<b>Species</b>	N/A
<b>Conjugate</b>	BSA
<b>Applications</b>	The antigen can be used in ELISA assay.
<b>Format</b>	Liquid
<b>Concentration</b>	Batch dependent - please inquire should you have specific requirements.
<b>Size</b>	1 mg
<b>Buffer</b>	0.01M pH7.4 PBS
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.