



Pyrimethanil [OVA] (DAG5639)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description Pyrimethanil is a member of the class of aminopyrimidines that is N-phenylpyrimidin-2-amine carrying two additional methyl substituents at positions 4 and 6. Pyrimethanil is a broad spectrum fungicide often applied to seeds. It inhibits methionine biosynthesis, thus affecting protein formation and subsequent cell division. Pyrimethanil works best on young fungus infestations.

Species N/A

Conjugate OVA

Applications The antigen can be used in ELISA assay.

Format Liquid

Concentration Batch dependent - please inquire should you have specific requirements.

Size 1 mg

Buffer 0.01M pH7.4 PBS

Storage Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.