



Pyrimethanil [KLH] (DAG5640)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description Pyrimethanil is a member of the class of aminopyrimidines that is N-phenylpyrimidin-2-amine carrying two additional methyl substituents at positions 4 and 6. Pyrimethanil is a broad spectrum fungicide often applied to seeds. It inhibits methionine biosynthesis, thus affecting protein formation and subsequent cell division. Pyrimethanil works best on young fungus infestations.

Species	N/A
Conjugate	KLH
Applications	The antigen can be used in ELISA assay.
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.