



Picoxystrobin [BSA] (DAG5644)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description Picoxystrobin is an enoate ester that is the methyl ester of (2E)-3-methoxy-2-[2-({[6-(trifluoromethyl)pyridin-2-yl]oxy}methyl)phenyl]prop-2-enoic acid. It is a cereal fungicide used to control a wide range of diseases including brown rust, tan spot, powdery mildew and net blotch. It is an aromatic ether, an enoate ester, an enol ether, an organofluorine compound, a member of pyridines and a methoxyacrylate strobilurin antifungal agent, has a role as a mitochondrial cytochrome-bc1 complex inhibitor and an antifungal agrochemical.

Species	N/A
Conjugate	BSA
Applications	The antigen can be used in ELISA assay.
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.