



Carazolol [BSA] (DAG5647)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Carazolol is a high affinity antagonist/partial inverse agonist (also referred to as a beta blocker) of the β -adrenergic receptor
Species	N/A
Conjugate	BSA
Applications	The antigen can be used in ELISA assay.
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.