



Metformin [BSA] (DAG5653)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Metformin is an agent belonging to the biguanide class of antidiabetics with antihyperglycemic activity. It is mainly used for the treatment of type 2 diabetes and can be used alone or in combination with sulfonylureas, thiazolidinediones or other hypoglycemic agents. Up to now, metformin is the most widely used medication for diabetes taken by mouth.
Species	N/A
Conjugate	BSA
Applications	The antigen can be used in ELISA assay.
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.