



Diclofenac [OVA] (DAG5654)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Diclofenac (marketed under many brand names, see below: Trade names) is a non-steroidal anti-inflammatory drug (NSAID) taken to reduce inflammation and as an analgesic reducing pain in certain conditions. The name is derived from its chemical name: 2-(2,6-dichloranilino) phenylacetic acid. In the United Kingdom, India, Brazil and the United States, it may be supplied as either the sodium or potassium salt, in China most often as the sodium salt, while in some other countries only as the potassium salt. Diclofenac is available as a generic drug in a number of formulations. Over-the-counter (OTC) use is approved in some countries for minor aches and pains and fever associated with common infections.
Species	N/A
Conjugate	OVA
Applications	The antigen can be used in ELISA assay.
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.