



Fluorescein [KLH] (DAG5659)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Fluorescein is an organic compound and can emit an intense green fluorescence in the presence of alkaline solutions, and is used forensically to detect traces of blood, in analytical chemistry as an indicator in silver nitrate titrations and in microscopy. Fluorescein has an absorption maximum at 494 nm and emission maximum of 512 nm (in water). The major derivatives are fluorescein isothiocyanate (FITC).
Species	N/A
Conjugate	KLH
Applications	The antigen can be used in ELISA assay.
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.