



Digoxin [KLH] (DAG5661)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Digoxin is a cardiac glycoside, which can inhibits the sodium potassium adenosine triphosphatase (ATPase) pump, thereby increasing intracellular calcium and enhancing cardiac contractility. It is named Lanoxin as a medication used to treat various heart condition, such as atrial fibrillation, atrial flutter, and heart failure. Digoxin is first extracted from the foxglove plant, <i>Digitalis lanata</i> , and now is the most effective and safe medicines needed in a health system.
Species	N/A
Conjugate	KLH
Applications	The antigen can be used in ELISA assay.
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.