



Human Cholinergic Receptor, Muscarinic 2 [His] (DAG625)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Human
Purity	>75% as determined by SDS-PAGE
Conjugate	His
Applications	WB, ELISA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	50 µg
Buffer	Neutral to slightly alkaline pH.
Preservative	None
Storage	2-8°C short term, -20°C long term

BACKGROUND

Introduction	The muscarinic acetylcholine receptor M2, also known as the cholinergic receptor, muscarinic 2, is a muscarinic acetylcholine receptor.
Keywords	Muscarinic acetylcholine receptor M2; Cholinergic receptor, muscarinic 2; CHRM2