



Human Chromodomain Helicase DNA Binding Protein 4 [His] (DAG626)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Human
Purity	>80% as determined by SDS-PAGE
Conjugate	His
Applications	WB, ELISA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	50 µg
Buffer	Neutral to slightly alkaline pH.
Preservative	None
Storage	2-8°C short term, -20°C long term

BACKGROUND

Introduction	Chromodomain-helicase-DNA-binding protein 4 is an enzyme that in humans is encoded by the CHD4 gene
Keywords	CHD4; Chromodomain helicase DNA binding protein 4