



6-MAM [BSA] (DAGA-034)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	BSA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	Distilled Water
Preservative	None
Storage	Short Term (≥ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.

BACKGROUND

Keywords	6-Monoacetylmorphine; 6-MAM
----------	-----------------------------