



# Recombinant Helicobacter pylori urease large subunit Protein (DAGA-061)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Antigen Description	Protein produced in genetically engineered E. coli cells. It covers Met 1 to LYS559.
Purity	> 90%
Conjugate	Unconjugated
Applications	N/A
Molecular Weight	63 KDa
Stability	At least 1 year at -20°C.
Format	Solution
Size	50 µg
Buffer	Solution at 1.0 mg/ml in 50 mM Tris-Cl, 0.1% SDS/1 mM EDTA (pH 7.5)
Preservative	None
Storage	Store at -20°C.

## BACKGROUND

Keywords	Helicobacter pylori urease large subunit antigen; H.pylori urease large subunit antigen
----------	---