



## 3-Methyl-Quinoxaline-2-Carboxylic acid [KLH] (DAGA-061K)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

|               |                                                                                       |
|---------------|---------------------------------------------------------------------------------------|
| Conjugate     | KLH                                                                                   |
| Applications  | Immunogen                                                                             |
| Format        | Liquid                                                                                |
| Concentration | Batch dependent - please inquire should you have specific requirements.               |
| Size          | 1 mg                                                                                  |
| Buffer        | 0.01M pH7.4 PBS                                                                       |
| Preservative  | None                                                                                  |
| Storage       | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. |

### BACKGROUND

|          |                                                                              |
|----------|------------------------------------------------------------------------------|
| Keywords | 3-Methylquinoxaline-2-carboxylic Acid;3-Methyl-Quinoxaline-2-Carboxylic acid |
|----------|------------------------------------------------------------------------------|