



DENV type 1, 4 (55 kDa) (DAGA-079)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Caused by one of four closely related virus serotypes of the genus Flavivirus, family Flaviviridae, each serotype is sufficiently different that there is no cross-protection and epidemics caused by multiple serotypes (hyperendemicity) can occur. In cell culture experiments and mice Morpholino antisense oligos have shown specific activity against Dengue virus. The E.coli derived recombinant 55kDa protein is genetically engineered peptide which is derived from Dengue Type-1 and 4 to be expressed as a fused envelope, each part in this fusion contains 170 a.a (positions 46-217), it is used in ELISA assay. This fusion protein is connected to a 6xHis Tag. Dengue Type-1 and 4 is purified by proprietary chromatographic technique.
----------------------------	---

Purity	Protein is >95% pure as determined by 12% PAGE (coomassie staining).
Conjugate	Unconjugated
Applications	ELISA
Molecular Weight	55 kDa
Stability	Stable at 4°C for 1 week. Please prevent freeze thaw cycles.
Format	Solution
Size	100 µg, 500 µg, 1 mg
Buffer	Phosphate buffered saline, pH-7.4.
Preservative	None
Storage	Stored below -18°C.
Ship	Ice Packs

BACKGROUND

Keywords

Dengue Virus 2+4; Dengue Virus Type 2, 4; Group IV; Flaviviridae; Flavivirus; Dengue virus; Dengue Virus; DENV; DENV type 1, 4
