



DENV type 2 Envelope protein (C-terminal, 22 kDa) (DAGA-099)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description Caused by one of four closely related virus serotypes of the genus *Flavivirus*, family *Flaviviridae*, each serotype is sufficiently different that there is no cross-protection and epidemics caused by multiple serotypes (hyperendemicity) can occur. In cell culture experiments and mice Morpholino antisense oligos have shown specific activity against Dengue virus. The E.coli derived recombinant 15kDa protein is genetically engineered peptide which is derived from C terminus of Dengue Type-2 envelope. The expressed protein contains neutralizing epitopes and receptor binding domain, it has been used in the study for vaccine development. A 6 x His Tag is fused at its C-terminus.

Purity Protein is >95% pure as determined by 12% PAGE (coomassie staining).

Conjugate Unconjugated

Applications N/A

Molecular Weight 15 kDa

Stability Stable at 4°C for 1 week. Please prevent freeze thaw cycles.

Format Solution

Size 100 µg, 500 µg, 1 mg

Buffer Phosphate buffered saline pH 7.4.

Preservative None

Storage Stored below -18°C.

Ship Ice Packs

BACKGROUND

Keywords

Dengue virus 2; Flavivirus; DENV Envelope Protein; DV 2 Envelope Protein; Dengue Virus 2 Envelope Protein; DENV type 2 Envelope protein
