



Recombinant DENV type 2 envelope protein (D2-D80-100, Fragment) [His] (DAGA-1007)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Recombinant 6xHis tagged DENV type 2 envelope protein (D2-D80-100, Fragment) (UniProtKB No. Q88629) was expressed.
Purity	>85%
Conjugate	His
Applications	ELISA, Immunogen
Recommended Usage	This item requires custom production and lead time is between 4-8 weeks. We can custom produce based on your specifications.
Format	Lyophilized
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg, 5 mg
Buffer	Each vial contains 100 µg of purified protein in PBS
Preservative	None
Storage	Store at 2-8 °C short term and -20°C long term. It is recommended to aliquot needed volumes and store at -20°C. Avoid repeated freeze/thaw cycles.

BACKGROUND

Keywords	DENV type 2; Dengue Virus type 2; Dengue Virus 1+2+3+4; Dengue Virus Type 1, 2, 3, 4; Group IV; Flaviviridae; Flavivirus; Dengue virus; Dengue Virus; DENV
-----------------	--

GENE INFORMATION

UniProt ID	Q88629
------------	--------
