



## Benzophenone [OVA] (DAGA-1000)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

Conjugate	OVA
Applications	ELISA, LFIA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	1 mg
Buffer	10mM PBS, pH7.4
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
Ship	Wet ice

### BACKGROUND

Introduction	Benzophenone is the organic compound with the formula (C <sub>6</sub> H <sub>5</sub> ) <sub>2</sub> CO, generally abbreviated Ph <sub>2</sub> CO. It is a white solid that is soluble in organic solvents. Benzophenone is a widely used building block in organic chemistry, being the parent diarylketone.
--------------	--