



Recombinant YFV Envelope protein (Fragment) [His] (DAGA-1074)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Recombinant 6xHis tagged YFV Envelope protein (Fragment) (UniProtKB No. Q89294) was expressed.
Purity	>85%
Conjugate	His
Applications	ELISA, Immunogen
Recommended Usage	This item requires custom production and lead time is between 4-8 weeks. We can custom produce based on your specifications.
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	50 mM Tris-HCl, 8M urea
Preservative	None
Storage	Store at 2-8 °C short term and -20°C long term. It is recommended to aliquot needed volumes and store at -20°C. Avoid repeated freeze/thaw cycles.

BACKGROUND

Keywords Yellow fever virus; Envelope protein; E protein

GENE INFORMATION

UniProt ID

Q89294
