



Recombinant YFV Non-structural protein 4A (Fragment) [His] (DAGA-1077)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Recombinant 6xHis tagged YFV Non-structural protein 4A (Fragment) (UniProtKB No. Q89303) was expressed.
Purity	>85%
Conjugate	His
Applications	ELISA, Immunogen
Recommended Usage	This item requires custom production and lead time is between 4-8 weeks. We can custom produce based on your specifications.
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg, 5 mg
Buffer	PBS, pH7.4
Preservative	None
Storage	Store at 2-8 °C short term and -20°C long term. It is recommended to aliquot needed volumes and store at -20°C. Avoid repeated freeze/thaw cycles.

BACKGROUND

Keywords Yellow fever virus; NS4A; Nonstructural protein

GENE INFORMATION

UniProt ID	Q89303
------------	--------
